# MediFlowTec™, MEDIFLOW® II, MEDIFLOW +

FLOW-METERING DEVICES  
(FLOW SELECTOR); INSTRUCTION FOR USE

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EN

## 1. FOREWORD

GCE Flow selectors are medical devices classified as class IIa according to the Medical Device Directive 93/42/EEC.

The Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon 15002 standard.

## 2. INTENDED USE

Mediflow are devices intended for the administration of the following medical gases in the treatment and care of patients:

- oxygen O<sub>2</sub>
- medical air AIR
- nitrous oxide N<sub>2</sub>O

The product is not intended to be use with air or nitrogen for driving surgical tool.

This flow selector is intended to be fitted to medical gas pipeline system terminal units in hospital or ambulance car or to quick connector outlet of medical regulator, with nominal pressures up to 5 bar.

**⚠ Using variants with maximum flow 50 l/min is necessary to be more careful and outlet of the flow selector MUST NOT be connected directly to the patient, for example via cannula.**

**⚠ The set value on the flow selector does NOT indicate real the actual flow of medical gas – the difference can be caused by different (missing) inlet pressure or any fault on outlet hose (blocked, fold...)**

## 3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

**⚠ Keep away the product and its associated equipment from:**

- heat sources (fire, cigarettes, ...),
- flammable materials,
- oil or grease
- water
- dust

**⚠ The product and its associated equipment must be prevented from falling over.**

 Use only the product and its associated equipment in well vented area.

 Always maintain oxygen cleanliness standards.

Before initial use the product should be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdraw from operation (for transport, storage). Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

EN

OPERATING CONDITIONS		STORAGE AND TRANSPORT CONDITIONS	
	-20/+60 °C		-30/+60 °C
	10/100%		10/100%
	600/1200 mbar		600/1200 mbar

 In case of storage at a temperature below -20 °C, do not operate the flow selector until it has been allowed to increase its temperature to a minimum of -20 °C.

## 4. PERSONNEL INSTRUCTIONS


The Medical Devices Directive 93/42/EEC states that product provider must ensure that all personnel handling the product are provided with the operating instructions & performance data.

 Do not use the product without properly familiarization of the product and its safe operation as defined in this Instruction for use. Ensure user is aware of particular information and knowledge required for the gas in use.

## 5. PRODUCT DESCRIPTION

Gas from the medical gas pipeline system is entering through the inlet connector (A). There are two basic versions: 1 – inlet connector is directly mounted on the flow selector, 2 – connection to the rail, including hose and inlet connector. Flow is then adjusted to the requested value by turning handwheel (B). The adjusted flow is visible in the window of the handwheel / on the pointer (C). Gas is delivered to the user through the flow outlet user connection (D).

EN

### A. INLET QUICK CONNECTOR

The flow selector is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe.

### B. HANDWHEEL

Requested gas flow can be adjusted by turning the handwheel to set the required value in the viewing window / pointer (calibrated in l/min). By the handwheel is possible to close the gas flow – adjusting value “0”.

### C. READING WINDOW

Adjusted value of the gas flow is visible in reading window / on the top side of the pointer.

### D. USER OUTLET

The “flow outlet” user connection is usually with connection for hose (hose nipple) or threaded (for humidifier).

This outlet is for supplying a gas flow (l/min) atmospheric pressure directly to a patient, e.g. through a cannula or a facemask.

 **The flow outlet shall not be used for driving any medical device.**

Note also that products colour (especially handwheel) might not follow any gas colour coding.

## 6. OPERATIONS

### 6.1. BEFORE USE

#### 6.1.1. VISUAL INSPECTION BEFORE FITTING

- Check if there is visible external damage of the flow selector (including product labels and marking). If it shows signs of external damages, remove it from service and identify its status.
- Visually check if the product is contaminated; if needed, for the flow selector, use the cleaning procedure detailed in this instruction for use.
- Ensure the handwheel is set on “0” position (turning anticlockwise).
- Ensure that connected hose nipple or humidifier to the outlet has correct dimensions.

 **The flow selector is dedicated only for use with the gas specified on its labelling. Never try to use for another gas.**

#### 6.1.2. FITTING TO MEDICAL GAS PIPELINE SYSTEM TERMINAL UNITS

The flow selector shall be connected in accordance with suitable type of terminal unit or with a regulator with an appropriate quick connector, see appendix nr.2.

#### 6.1.3. FUNCTIONAL TEST

- Check that there is a gas flow at each setting higher than “0” set by the handwheel (for instance, by listening of the sound of gas flow or checking presence of bubbles in a humidifier).
- Set the handwheel back to the close position – value “0”.

 **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

### 6.2. CONNECTION & USE OF OUTLET

#### 6.2.1. LIST OF RECOGNIZED ACCESSORIES

To be connected to the flow outlet:

- Humidifier, breathing mask or cannula.

 **Before connecting any accessory or medical device to the flow selector, always check that it is fully compatible with the product connection features & the product performances and there is no negative influence to outlet parameters of the flow selector.**

#### 6.2.2. FLOW OUTLET CONNECTION

 **When connecting any accessory to the flow outlet make sure that it is not connected to the patient before operating the product.**

- Ensure the hose/humidifier is compatible with the flow outlet feature.
- Push the hose onto the flow selector flow outlet/screw onto the humidifier.
- Ensure the hose/humidifier is well engages.

### 6.2.3. FLOW SETTING

- Connect the flow selector to the source of medical gas.
- Set requested medical gas flow by the handwheel.

 Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position or in zero position.

 Always ensure that the flow control knob has engaged and not placed between two settings otherwise the flow selector will not deliver the correct flow of medical gas.

 During use ensure that connected hose is not mechanically damaged, throttled or broken, which would reduce or break supply of gas to the patient.

 The medical gas flow rate must be prescribed and administered by a clinically trained user.

 Fluctuating inlet pressure and ambient temperature can have influence on accuracy of flow measurement.

### 6.3. AFTER USE

- Before disconnecting flow selector from the medical pipe line system terminal unit or medical regulator, set the handwheel on close position - value "0".
- Remove connected accessories from user outlet.
- Disconnect flow selector from user outlet.

## 7. CLEANING

Remove general contamination with a soft cloth dampened in oil free oxygen compatible soap water & rinse with clean water.

Disinfections can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and that they are compatible with the product materials (including labels) and gas.

 Do not use cleaning solutions containing ammonia!

 Do not immerse in water or any liquid.

 Do not expose to high temperature (such as autoclave).

 To apply the cleaning solution, do not spray it, as the spray may enter the inner parts of product and cause contamination or damage.

 Do not use pressure water as it could damage or contaminate the product.

 Do not use the product under any circumstances of the inner parts of the product have been contaminated. It must be withdrawn from service.

## 8. MAINTENANCE

### 8.1. PRODUCT LIFE TIME

Depending on the date of production you can find on the product the marking in one of two versions are easy to see and easy to differentiate.

#### 8.1.1. MARKING PRIOR TO IMPLEMENTING UDI CODING

Form of nine digits serial number stamped on the product is following:  
YYMMXXXXX

YY: year of production

MM: month of production

XXXXX: sequence number

Example: serial number 190300521 shows the product produced in March 2019, with sequence number 521.

An additional 2D code (for internal usage in production process in GCE) is also part of product's marking .

#### 8.1.2. MARKING AFTER UDI CODING IMPLEMENTATION

There are 3 marked identifiers on the medical device surface:

14 numbers with prefix (01): GTIN - Global Trade Identification Number (identification of the selector type)

9 numbers with prefix (21): serial number in format YYMMXXXXX

6 numbers with prefix (17): the end of the product's life time in format YYMMDD

Example: marking number (17) 290227 shows the medical device with end of lite time on 27th of February 2029.

As the integral part of UDI marking is a 2D code, which consists of the same data. The batch number is included in the 2D code (prefixed (10) ) and it is for manufacturer use.

#### 8.1.3. MAINTENANCE

No special maintenance or service, apart from the tests before use, is needed. However, to make sure the product is in good working order, it would be good if the owner/cylinder distributor performs the checks (see 6.1) himself on a regular basis (ex. once every second year) and/or at every cylinder exchange. This is just for the owner to ensure that the product works well, especially in cases where the user has some health problems and is not able to check the product himself properly.

#### 8.1.4. PRODUCT LIFE TIME AND WASTE MANAGEMENT

Maximum life time of the product is 10 years. At the end of the product's life time, the product must be withdrawn from service (max.10 years).

The owner of the device shall prevent the reuse of the product and handle the product in compliance with “Directive of European Parliament and Council 2008/98/EC on waste“. In accordance to Article 33 of REACH GCE, s.r.o. as responsible manufacturer shall inform all customers if materials containing 0.1% or more of substances included in the list of Substance of Very High Concern (SVHC). The most commonly used brass alloys used for bodies and other brass components contain 2-3% of lead (Pb), EC no. 231-468-6, CAS no. 7439-92-1. The lead will not be released to the gas or surrounding environment during normal use.

After end of life the product shall be scrapped by an authorized metal recycler to ensure efficient material handling with minimal impact to environment and health. To date we have no information that indicates that other materials containing SVHC of concentrations exceeding 0.1% are included in any GCE product.

## **8.2. SERVICE, REPAIRS AND MAINTENANCE**

### **8.2.1. SERVICE**

GCE recommends that a product periodic inspection is undertaken at least once a year to do leakage test and check proper functionality of the flow selector.

### **8.2.2. LEAKAGE TEST**

- After connection to medical gas pipeline system terminal unit or to medical regulator, visually check possible leakages:

#### **LEAKAGE TEST OF THE FLOW SELECTOR FOR FLOW ADJUSTMENT:**

- 1.** Adjust by handwheel the value „0“. When closing don't use excessive torque.
- 2.** Connect the flow outlet and hose nipple with a hose, the end of the hose immerse into a jar filled with water.
- 3.** Connect the flow selector to the source of compressed air max. 5 Bar.
- 4.** The test: leave the flow selector 1min under pressure. Number of bubbles shall not be higher than 3 bubbles per 1 min.
- 5.** When test ready, disconnect the flow selector from the source of compressed air.

#### **EXTERNAL LEAKAGE TEST:**

- 1.** Adjust by handwheel of the flow selector the value „0“. When closing don't use excessive torque.
- 2.** Close the outlet connection of the flow selector using a plug, the plug has to have the same thread as the outlet connector of the flow selector or, in the case the hose nipple used, use the hose with blind plug.
- 3.** Connect the flow selector to the source of compressed air max. 5 Bar.
- 4.** Adjust by handwheel of the flow selector the maximal value. When opening don't use excessive torque.
- 5.** Possible leakage you will detect by application of soap water on product joins.

6. After test disconnect the flow selector from the source of compressed air, remove the plug from outlet connection and carefully clean the surface of the flow selector, see chapter 7 (cleaning).

**⚠** If any leakage is detected, use the procedure in chapter 6.3 and return the pro-duct to GCE for service.

**⚠** The product shall not be disassembled or repaired when connected to the source of the pressure.

EN

#### REPAIRS ACTIVITIES COVER THE REPLACEMENT OF THE FOLLOWING DAMAGED OR MISSING COMPONENTS:

- quick connector (probe),
- flow head (valid only for MFII and MF+)
- plastic cover and handwheel (valid only for MFT)

**⚠** The repairs shall be carried out by a GCE authorised person only.

Any product sent back to a GCE authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair, maintenance). For product to be repaired a short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. The following parts can be replaced only:

- Caps,
- Stickers,
- Hose nipple (including o-ring),

**⚠** To use suitable component, contact our service organization.

**⚠** All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.

**⚠** All seals and o-rings must be kept in dry, dark and clean environment by the owner and the user during the entire product life time.

**⚠** Use only original GCE components!

## 9. GLOSSARY

	Consult instruction for use		Suitable for Home care use
	Caution		Suitable for Hospital care use
	Keep away from heat and flammable material		Suitable for Emergency care use
	Keep away from oil and grease		Serial number

	Humidity limit	<b>REF</b>	Catalogue number
	Temperature limit	<b>LOT</b>	Batch code
	Keep dry		Fragile, handle with care
	Date of manufacture		Manufacturer
	Use by date	<b>P<sub>1</sub></b>	Inlet pressure
	Inlet parameter		Outlet parameter

## 10. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

### APPENDIX:

Nr 2 - Quick coupling feature and connecting/disconnecting procedure

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Czech Republic              © GCE, s.r.o.