# VARIMED

---

**HIGH PRESSURE REGULATORS**

EN

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ES

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PT

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DE

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NL

**DETENDEURS HAUTE PRESSION**

FR

INSTRUCTION FOR USE

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MODE D'EMPLOI

## 1. FOREWORD

GCE Medical Regulators are medical devices classified as class IIb according to the Medical Device Directive 93/42/EEC.

Their Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon ISO 10524-2 standard.

## 2. INTENDED USE

Pressure regulators equipped with gas cylinder inlet connection are to be used in emergency as back up gas supply for the hospital pipeline system. **They are NOT designed for direct use with patient.** They are intended for the administration of the following medical gases:

- oxygen;
- nitrous oxide;
- air for breathing;
- helium;
- carbon dioxide;
- xenon;
- specified mixtures of the gases listed.
- air or nitrogen to power surgical tools.

## 3. OPERATIONAL, TRANSPORT AND STORAGE SAFETY REQUIREMENTS

**⚠ Keep the product and its associated equipment away from:**

- All sources of heat
- Flammable materials
- Oil or grease (including all hand creams)
- Water
- Dust.

**⚠ The product and its associated equipment must be prevented from falling over.**

**⚠ Always maintain oxygen cleanliness standards.**

**⚠ Use only the product and its associated equipment in a well ventilated area.**

Before initial use the product should be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdrawn from operation (for transport, storage).

Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

## 4. PERSONNEL INSTRUCTIONS AND TRAINING


The Medical Devices Directive 93/42/EEC states that the product provider must ensure that all personnel using the product are provided with the instructions for use and are fully trained in the use of the equipment.

**⚠ Do not use the product without proper training! Trainees must be trained by an experienced person who has been authorised by the manufacturer and has an appropriate education, knowledge and experience.**

For further information about training of responsible personnel, please contact GCE.

## 5. PRODUCT DESCRIPTION

Pressure regulator reduces gas pressure from the cylinder or in the distribution system. Gas passes the cylinder through the pressure regulator to the gas outlets.

EN

### **A - Inlet stem**

Pressure regulator is fitted to the cylinder shut-off valve or to the distribution system by means of an inlet stem. The stem can have a nut with female or male thread. The inlet stem includes a filter.

### **B - Inlet pressure indicator**

The regulator is fitted with a pressure indicator which is intended for cylinder gas content indication only, not for measuring purposes.

### **C - Outlet pressure indicator**

Outlet pressure indicator indicates outlet pressure of the valve.

### **D - Pressure outlet**

Pressure regulator is fitted with a pressure outlet. The pressure outlet is the outlet directly from the low-pressure chamber. Two types of the pressure outlet can be used:

Pressure outlet I – is fitted with a gas specific medical quick connector also called „quick coupler“. The quick connector self seals when the male probe is disconnected.

Pressure outlet II – is fitted with a threaded connector.

### **E - Control knob**

Control knob is used for adjustment of outlet pressure (if applicable).

### **Pressure relief valve**

Pressure relief valve serves for reduction of overpressure in the regulator. It is fitted on the rear part of the regulator body.

## **6. OPERATIONS**

### **6.1. BEFORE USE**

#### **6.1.1. Visual inspection before use**

- Check the product (including product labels and marking) and the gas cylinder for damage. If damaged, remove it from service and identify its status suitably.
- Visually check if the product or the medical gas cylinder is contaminated; if needed, for the regulator, use the cleaning procedure detailed in this document (if required for the cylinder, refer to the gas cylinder manufacturer cleaning recommendation).
- Check if the product service is due or that the total life time of the GCE product has not been exceeded according to the serial number stamped on the rear part of the regulator (see section 9). If service or life time has been exceeded, remove the regulator from service and suitably identify its status.
- Ensure that the product inlet stem is compatible with the medical cylinder valve (gas/thread type).
- Check the presence & the integrity of inlet stem seals / correct size of seal.

**⚠ The product is intended only for use with the gas specified on the product. Never use it for other gas.**

#### **6.1.2. Connection to medical cylinder valve**

- Ensure that you have the correct medical gas cylinder. Refer to your gas supplier for medical gas data sheets.
- Remove protective cap from cylinder valve, direct the cylinder valve away from personnel or naked flames and briefly open the valve to ensure that water or foreign matter is ejected from the valve outlet.
- Before fitting to the cylinder ensure that the valve outlet and the regulator inlet are clean and free from oil and grease.
- Fit the regulator into the cylinder valve and secure in position.
- Cylinders should be located where they will not be contaminated with oil or grease and securely supported in a vertical position.

**⚠ Fitting the regulator with too high tightening torque to the cylinder valve may result in damage.**

**⚠ During fitting to the cylinder valve, do not apply torque/load to any other parts of the product.**

#### **6.1.3. Leakage check before use**

- Open the cylinder valve slowly by turning the hand wheel in anticlockwise direction approx. 1 to 1.5 turns.

**⚠ Sudden opening of the cylinder valve could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the cylinder valve could reduce actual flow delivered.**

- Visually and by listening check possible leakages on:
  - regulator inlet connection to cylinder valve,
  - pressure indicator to regulator body,

- pressure relief valve vent holes.
- Turn off the cylinder valve by turning the hand wheel in clockwise direction to „stop“ position. Do not use excessive force

**⚠ If any leakage is detected, use the procedure in chapter 6.3 and return the product for service.**

#### **6.1.4. Functional test before use**

- Open the cylinder valve – position “ON”.
- Check that the pressure gauge indicates pressure.
- Turn off the cylinder shut-off valve. Do not use excessive force.
- To test the regulator functionality, check gas flow through the regulator. For regulators with a quick coupler outlet, connect and disconnect the male QC probe. For regulators with threaded outlet, gas flows already in opening the cylinder valve.

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**⚠ Sudden opening of the cylinder valve could result in a danger of fire or explosion arising from oxygen pressure shocks. Insufficient opening of the cylinder valve could reduce actual flow delivered.**

### **6.3. AFTER USE**

- Turn off the cylinder valve by turning the hand wheel in clockwise direction to “stop” position. Do not use excessive force.
- Vent the regulator.
- Ensure the pressure indicator does not show any residual pressure.
- Disconnect all connected equipment from the user outlets.

**⚠ Never leave regulators under pressure unless they are in use!**

## **7. CLEANING**

Remove general contamination with a soft cloth damped in oil free oxygen compatible soapy water & rinse with clean water.

Disinfection can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and that they are compatible the product materials (including labels) and gas.

**⚠ Do not use cleaning solutions containing ammonia!**

**⚠ Do not immerse in water or any liquid.**

**⚠ Do not expose to high temperature (such as autoclave).**

## 8. MAINTENANCE

### 8.1. SERVICE

1. Regulators should be regularly inspected for signs of mechanical damage or leakage. If incorrect function or leakage is suspected, the regulator should be immediately returned to the manufacture for service exchange or replacement.
2. Pressure regulators are delicate precision instruments and no attempt should be made to remove or exchange components from them.
3. Washers should be carefully examined to ensure they are undamaged and have not gone brittle or sticky. Only original manufactures' washers must be used. Regulators should be removed from service and fully leak tested annually, and returned for service exchange after a maximum of six years service. When used in rugged and severe conditions these times should be reduced.
4. Oil or grease must never be allowed in contact with the regulator, or cylinder valve, or with any accessories used with the equipment.
5. PTFE tape must not be used on this regulator.

### 8.2. LIFE TIME

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be withdrawn from service. The owner shall put in place a relevant procedure to ensure the product cannot be used again.

### 8.3. REPAIRS

Repairs activities cover the replacement of the following damaged or missing components:

- Inlet stem,
- Indicator,
- Pressure relief valve,
- Quick coupler/Outlet stem

The repairs shall be carried out by a GCE authorised person only.

Any product sent back to a GCE authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair, overall maintenance). For product to be repaired a short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. The following parts can be replaced only:

- Caps,
- Hose nipple (including o-ring)
- O-ring for inlet stem

**⚠ All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

**⚠ All seals and o-rings must be kept in dry, dark and clean environment by the owner and the user during the entire product life time.**

**⚠ Use only original GCE components!**

## 9. GLOSSARY

	Consult instruction for use		Suitable for Hospital care use
	Caution		Keep away from oil and grease
	Keep away from heat and flammable material		Upper and lower humidity limit
	Upper and lower temperature limit		Product serial number YY MM XXXXX Y: manufacturing year M: manufacturing month X: production number
	Keep dry		
	Fragile		Reference number
	Use by date		Batch number
	Inlet pressure range		GCE number
	Date of manufacture		Manufacturer
	Outlet parameter		Inlet parameter

EN

## 10. WARRANTY

The Standard Warranty period is two years from date of receipt by the GCE Customer (or if this is not known 2 years from time of the product manufacture shown on the product).

The standard warranty is only valid for products handled according to Instruction for use (IFU) and general industry good practice and standards.

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
 2460

GCE Group is one of the world's leading companies in the field of gas control equipment. The headquarters are in Malmö, Sweden, and the two major supply units are located in Europe and Asia.

The company operates 15 subsidiaries around the world and employs more than 850 people. GCE Group includes four business areas –Cutting & Welding, Process Applications, Medical and High Purity. Today's product portfolio corresponds to a large variety of applications, from single pressure regulators and blowpipes for cutting and welding to sophisticated gas supply systems for medical and electronics industry applications.

**Gas Control Equipment**

GCE world-wide: <http://www.gcegroup.com>

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Czech Republic            © GCE s.r.o.