



NOTE :
 Used without flap when packing from BOM1
 Used with flap when packing from BOM6

1	1	Hot dipped galvanized cradle	1004361	16001572	Fe 360
POS:	Qty.	Description	Material No:		
If the product is requested to be manufactured differently, this MUST be approved by VIKING LIFE-SAVING EQUIPMENT A/S and a specification of change must be filled in if the alteration is approved. Unless otherwise specified the product must always be clearly labelled with Viking part no., manufacturer ID and manufacturing year.					
REVISION NO:	REVISION DESCRIPTION:			REVISION DATE:	APPROVAL INITIALS:
02	Add service hatch on container for BOM6			01.04.19	PSOM
MATERIAL NO:	STATUS:	UNITS:	CHECK DATE/INITIALS:		
1042411	-	mm	27.05.16/ EB		
BASIC MATERIAL:	SURFACE TREATMENT:		SHEETSIZE:	SCALE:	SHEET NO:
-	-		A3	-	-
VIKING LIFE-SAVING EQUIPMENT A/S Sædding Ringvej 13, DK-6710 Esbjerg V Tel.: +45 7611 8100 Fax.: +45 7611 8101 E-mail: viking@viking-life.com Web: www.VIKING-life.com			PROJECTION TYPE:		DESCRIPTION:
			 FIRST ANGLE PROJECTION GENERAL TOLERANCE: -		
Copyright of VIKING LIFE-SAVING EQUIPMENT DENMARK whose property this document remains. No part thereof may be disclosed, copied, duplicated or in any other way made use of except with the approval of VIKING LIFE-SAVING EQUIPMENT DENMARK			DRAWING NO.:		
			42002786		

Description of change